

**The Devil is in the Fitting Details**

Why all NAL (or DSL) targets are not created equal

Mona Dworsack-Dodge, Au.D.  
Senior Audiologist- Otometrics, DK  
August, 2012  
Audiology Online 20961

---

---

---

---

---

---

---

---

**What they share in common...**

- Evidence based prescriptions with a host of literature validating them as tools for best practice
- Use of speech-like signals as stimulus of choice
- Differing targets for soft, average, and loud inputs
- Provide a consistent fitting approach regardless of make, model, manufacturer of device

---

---

---

---

---

---

---

---

**NAL-NL2**

- **PREMISE: Loudness EQUALIZATION**
  - Aims to equalize perception of loudness over a range of frequencies instead of having low frequencies dominate loudness
  - maximizing predicted speech intelligibility while constraining loudness
- **EVOLUTION:**
  - From LINEAR approach to COMPRESSION approach.
  - From INSERTION-GAIN (tones) to include AIDED Gain (speech) targets
- **VALIDATION/EVIDENCE:** modifications from NAL-NL1 included patient preference/comfort findings

Johnson, E. and Dillon, H., 2011

---

---

---

---

---

---

---

---

### Desired Sensation Level Multistage Input/Output (DSL m[i/o])- DSL v5.0

- PREMISE: Loudness NORMALIZATION
  - Aims to restore, at each frequency loudness perception of the hearing impaired listener to that of a normal hearing listener
  - Goals of avoiding loudness discomfort, providing audibility of speech across a wide range of input levels, and accommodating the prescriptive targets for both quiet and noisy environments, as well as for infants versus children versus adults
- EVOLUTION:
  - \*\*From pediatric-focus (earlier versions) to both pediatric and adult versions.\*\*
  - Modifications for 'noise programs' made
  - Correlates with data on Preferred Listening Levels (PLL)

Johnson, E. and Dillon, H., 2011

---

---

---

---

---

---


---

---

---

---

### Where are they now? Soft (55dB) input




---

---

---

---

---

---


---

---

---

---

### Where are they now? Average (65) input




---

---

---

---

---

---


---

---

---

---

## Aided Response - Average Speech




---

---

---

---

---

---

---

---

---

---

## INTER-SYSTEM DEVIATIONS

---

---

---

---

---

---

---

---

---

---

### Whose NAL-NL fitting method are you using?

By Todd A. Ricketts and H. Gustav Mueller

There are many factors to consider in selecting hearing aid gain and output for a patient. Do you choose the fitting that maximizes audibility? The fitting that optimizes speech understanding? The fitting that "sounds the best" to the patient? The fitting recommended by the manufacturer?

Using these or other possible "gold

correction factors during the programming and verification process. With infants and young children, it may be necessary to measure the RECD and then use these values to derive estimated real-ear output. With adults, equipment and protocols are available to conduct real-ear verification measures directly.

**ABBREVIATIONS:**  
 DTE: Desired to use  
 DSE: Desired Sensation Level  
 SCL: Single-Pole Equalization Filter  
 NAL: NAL-NL2  
 NAL-NL1: National Acoustic Laboratories  
 NAL-NL2  
 FIDELITY  
 FIDELITY: Prescription of gain and output  
 RECD: Real-ear aided gain  
 RECD: Prescription of RECD

Ricketts, T.; Mueller H.G. (2009).

---

---

---

---

---

---

---

---

---

---

### Five factors which can affect accuracy of the verification process

- Target adjustments for the hearing aid/fitting type
- Conversions used to display hearing loss and targets in ear canal SPL
- Hearing aid-specific interactions with the measurement signal
- The level and shape of the input signal
- The analysis of the measured signal

Ricketts, T.; Mueller H.G. (2009).

---

---

---

---

---

---

---

---

### Inter-system deviations

- Compared NAL-NL1 targets for open and occluded fittings across 3 verification systems, using "system defaults."
- Randomly select one of the systems to fit to target and then measure (without adjustments) on the other two systems, and repeated cycle with the other two systems as "baseline."
- Calculated a deviation from NAL-NL1 targets for each system.
- Found results on Verifit and MedRX systems to be similar, but the Fonix system resulted in a desired fitting that was approximately 3-4 dB lower in the low frequencies, and up to 8-10 dB higher in the high frequencies.

Ricketts, T.; Mueller H.G. (2009).

---

---

---


---

---

---

---

---



**Figure 2.** The average intersystem deviation from NAL-NL1 targets across three probe-microphone systems for one closed fitted instrument. All values were calculated as described in Figure 1.

Ricketts, T.; Mueller H.G. (2009).

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

---

---

- ### Five factors which can affect accuracy of the verification process
- Target adjustments for the hearing aid/fitting type
  - Conversions used to display hearing loss and targets in ear canal SPL
  - Hearing aid-specific interactions with the measurement signal
  - The level and shape of the input signal
  - The analysis of the measured signal
- Ricketts, T.; Mueller H.G. (2009).

---

---

---

---

---

---

---

---

---

---

**NAL**

---

---

---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

## NAL NL2

- *The only mandatory inputs from the user are the air conduction hearing thresholds at 500 Hz and 2 kHz.*

---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

REAL-ML2

Client #  Name  DOB: 01-01-1990 Left Right Copy data to the other ear

Client Data    Autological Input    Selection Screen    Target Screen    Speech Screen

Hearing Aid Type

- BIC
- CIC
- CMC
- CICX

Number of Hearing Aids

- Unilateral
- Bilateral

Compression / Threshold

Fitting Depth

- Full
- Standard
- Deep

Acoustics

Tube

415

Sealing

Standard

Suggestions

Taking Selection i OK

Best Selection - This alert is likely to notice the occlusion effect with this aid.

Compression Speed

- Slow
- Fast
- Dual / Intermediate / Adaptive

Input quantity

- @45 dB
- @60 dB
- @75 dB in noise +20dB
- @85 dB in noise +30dB

Gain word thresholds: all HL

250Hz	500Hz	1kHz	2kHz	4kHz	8kHz	16kHz	32kHz
AC	15	15	40	45	55	60	65

DC

Use basic compilation values

Hearing Level

Hearing Threshold (dB HL) vs Frequency (Hz)

---

---

---

---

---

---

---

---

---

---

---

---

REAL-ML2

Client #  Name  DOB: 01-01-1990 Left Right Copy data to the other ear

Client Data    Autological Input    Selection Screen    Target Screen    Speech Screen

Hearing Aid Type

- BIC
- CIC
- CMC
- CICX

Number of Hearing Aids

- Unilateral
- Bilateral

Compression / Threshold

Fitting Depth

- Full
- Standard
- Deep

Acoustics

Tube

415

Sealing

Standard

Suggestions

Taking Selection i OK

Best Selection - This alert is likely to notice the occlusion effect with this aid.

Compression Speed

- Slow
- Fast
- Dual / Intermediate / Adaptive

Input quantity

- @45 dB
- @60 dB
- @75 dB in noise +20dB
- @85 dB in noise +30dB

Gain word thresholds: all HL

250Hz	500Hz	1kHz	2kHz	4kHz	8kHz	16kHz	32kHz
AC	15	15	40	45	55	60	65

DC

Use basic compilation values

Hearing Level

Hearing Threshold (dB HL) vs Frequency (Hz)

---

---

---

---

---

---

---

---

---

---

---

---

REAL-ML2

Client #  Name  DOB: 01-01-1990 Left Right Copy data to the other ear

Client Data    Autological Input    Selection Screen    Target Screen    Speech Screen

Hearing Aid Type

- BIC
- CIC
- CMC
- CICX

Number of Hearing Aids

- Unilateral
- Bilateral

Compression / Threshold

Fitting Depth

- Full
- Standard
- Deep

Acoustics

Tube

415

Sealing

Standard

Suggestions

Taking Selection i OK

Best Selection - This alert is likely to notice the occlusion effect with this aid.

Compression Speed

- Slow
- Fast
- Dual / Intermediate / Adaptive

Input quantity

- @45 dB
- @60 dB
- @75 dB in noise +20dB
- @85 dB in noise +30dB

Gain word thresholds: all HL

250Hz	500Hz	1kHz	2kHz	4kHz	8kHz	16kHz	32kHz
AC	15	15	40	45	55	60	65

DC

Use basic compilation values

Hearing Level

Hearing Threshold (dB HL) vs Frequency (Hz)

---

---

---

---

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

---

---

---

---

Why does venting not appear to influence Real Ear SPL targets?

---

---

---

---

---

---

---


---

---

---

---

---



Taken from the NAL-NL2 Algorithm implementation documentation

---

---

---

---

---

---

---

---


---

---

---

---





---

---

---

---

---

---

---

---

Why does venting not appear to influence Real Ear SPL targets?

*“There is no reason why changing the size of the venting should change the ear canal SPL that is optimal for a person at any frequency.”*

Dillon, H. 2006. What's new from the NAL in hearing aid prescriptions? *Hearing Journal* 10:10-16.

---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

---

---

---

---

Summary of variables and their effect

Category	Parameter	REHG	REAG	Ter CG
Client information	Air conduction thresholds	N	N	N
	Bone conduction thresholds	N	N	N
	Gender	N	N	N
	Age	N	N	N
	Hearing aid experience	N	N	N
	RECD	N	N	N
	RECU	N	N	N
	REDD	N	N	N
	REHG	N	(N)	N
	Language type	N	N	N
Hearing aid information	Bilateral - unilateral	N	N	N
	Compressions opened	N	N	N
	Limiting type	N	N	N
	Tubing	N	N	N
	Venting	N	N	N
	No of channels (a-b stimuli)	N	N	N
	No of channels (b-b stimuli)	N	N	N
Options	Wideband CT	N	N	N
	Limiting type	N	N	N
	Reference position	N	N	N
	Reference orientation	N	N	N
	Target type for coupler gain	N	N	N

Taken from the NAL-NL2 Algorithm implementation documentation

---

---

---

---

---

---

---

---

---

---

---

---

Ran: Audiometry Items

Ear	Left
Targets	DSL adult
Age	Adult
Transducer	Insert+Foam
Threshold	Entered
Bone conduction	N/A
UCL	Average
RECD	Average
Binaural	No
<input type="button" value="Cancel"/> <input type="button" value="Continue"/> <input type="button" value="Help"/>	

---

---

---

---

---

---

---

---


---

---

---

---

Rang Items



The screenshot shows a software interface for configuring hearing aid settings. It features several sections with frequency response curves and gain sliders. The sections include: 'AID for 300 Hz to 8 kHz', 'AID for Threshold', 'AID for 250 Hz', 'AID for 500 Hz', 'AID for 1000 Hz', and 'AID for 2000 Hz'. Each section has 'Left' and 'Right' ear settings. The curves show gain in dB across different frequencies. At the bottom, there are 'OK' and 'Cancel' buttons.

---

---

---

---

---

---


---

---

---

---

Range of options on PMM systems



The screenshot shows a 'Fitting prescription settings' dialog box. It is divided into 'Default prescription' and 'Advanced settings' sections. The 'Advanced settings' section includes: 'Transducer' (Headphone), 'Orientation for B110/RE110' (45 Deg), 'Hearing instrument style' (Behind the ear), 'Hearing aid type' (Linear), 'Tubing' (Elder\*), 'Earthing' (Off), 'Number of hearing aids' (Binaural), 'Threshold (dB)', 'Audogram type' (Use Test values), and 'Vented value' (Occluded). At the bottom are 'OK' and 'Cancel' buttons.

---

---

---

---

---

---


---

---

---

---

Range of options on PMM systems



The screenshot shows a 'Fitting Details' dialog box. It contains fields for: 'Target HLE' (NLS-NL2), 'Fitting Mode' (Coupler), 'Applied RELEG' (Predicted (NLS-NL2)), 'H.L. Type' (BTE), 'Verbing' (Occluded), 'Amplification' (Standard), 'Date of Birth' (01-01-1990), 'Gender' (Unknown), 'Applied RECD' (Predicted (NLS-NL2)), 'Transducer' (Headphone), 'Use Bone Conduction' (No), 'Applied HCL' (Predicted), 'Fabric' (AT3), 'Coupler' (PS1), 'Limiting' (Wideband), 'Coupler type' (2cc Coupler), 'No. of Channels' (1), 'Capacitor' (From Top), 'Compression' (52), 'RELEG Orientation' (0), 'Compression speed' (Fast), 'Language type' (Horizontal), 'Fitting Epoch' (Standard). At the bottom are 'Apply', 'Close', and 'Help' buttons.

---

---

---

---

---

---

---

---

---

---

**CONVERSION FACTORS**

---

---

---

---

---

---

---

---

**NAL-Potential sources of inter-system variance**

- Defaults assumed by each system may vary
- HL to SPL conversions (RECD, REUG, REDD) for each system may vary
- Transducer selection may not be aligned

---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

### HL to SPL Conversion in AURICAL

- In correspondence with the transducer selection, the RECD, REUG and any parameter impacting their calculation influences the HL/SPL conversion. 3 basic formulas apply (ANSI 3.6:2004):

---

---

---

---

---

---

---

---

### For HL to SPL (real ear)

- For inserts: Threshold [dB SPL] = Threshold [dB HL] + Applied RECD (mold or tip) + insert phone RETSPL (HA1 or HA2)
- For sound field: Threshold [dB SPL] = Threshold [dB HL] + Applied REUG (0, 45, or 90 degree) + Sound field RETSPL (0, 45, or 90 degree)
- For headphones: Threshold [dB SPL] = Threshold [dB HL] + real ear-to-6cc transform + RETSPL
- (For SPL to HL, the above formulas apply in reverse.)

---

---

---

---

---

---

---

---

### For SPL (real ear) to SPL (2cc)

- Threshold [dB SPL 2cc] = Threshold [dB SPL real ear] - Applied RECD

---

---

---

---

---

---

---

---

### SIGNALS AND TARGETS

---

---

---

---

---

---

---

---

### NAL-Potential sources of inter-system variance

- NAL doesn't provide REAR or coupler SPL targets
  - Derived by adding the AG to the input signal
  - Customized REAR/ Coupler SPL targets may differ amongst systems
    - E.G. Otosuite applies the known signal spectrum in order to achieve appropriately differing output targets for various signals

---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

### NAL-Potential sources of inter-system variance

- Bilateral/Unilateral settings may not be aligned
- Experience settings may not be aligned
- Proprietary Coupler Output for Flat Insertion Gain (CORFIG) Hawkins D, Cook J. Hearing aid software predictive gain values: how accurate are they? *Hear Jour.* 2003;56(7):26-34.

---

---

---

---

---

---

---

---

DSL

---

---

---

---


---

---

---

---

### The Desired Sensation Level (DSL) Method Version 5: $DSL^m[i/o]$



---

---

---

---

---

---

---

---



### DSL v5.0b m[i/o] -Required Inputs

- Age
- Client Type
- Transducer
- Audiogram including BC
- LDLs/UCLs
- Handling of ABR
- REUG and measurement conditions
- RECD and measurement conditions
- Linear vs. WDRC
- HI Style
- Monaural/Binaural
- # Channels
- Vent Size
- Signal type and multi-level
- Quiet vs Noise program

---

---

---

---

---

---

---

---

The screenshot shows the 'DSL m[i/o] 2.0D Application' window. It includes sections for 'Patient Description', 'Assessment Information', and 'Hearing Instrument Data'. The 'Assessment Information' section contains a table for Thresholds (HL) and a table for RECD. The 'Hearing Instrument Data' section includes fields for 'Number of Channels' and 'Crossover Frequency'.

Threshold (HL)	250	300	350	40	45	50	55	60	65	70
LDL (HL)										
UCL (HL)										
RECD										
REUG										

---

---

---

---

---

---

---

---

The screenshot shows the 'DSL m[i/o] 5.0B Application' window. It includes sections for 'Program 1' and 'Program 2'. Each program section contains a table for 'Verification Signal' and 'Output Format'. The 'Verification Signal' table includes columns for 'Level' and 'RECD'.

Level	50	55	60	65	70	75
Level 1: 52	50	55	60	65	70	75
Level 2: 60	60	65	70	75	80	85
Level 3: 74	74	79	84	89	94	99
RECD 90	90	95	100	105	110	115

---

---

---

---

---

---

---

---

What influence does venting have for DSL?

---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

What about HI Type?

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

---

---

---

---

### Range of options on PMM systems




---

---

---

---

---

---

---

---

---

---

### DSL-Potential sources of inter-system variance

- Defaults assumed by each system may vary
- HL to SPL conversions (RECD, REUG, REDD) for each system may vary
- UCL applied by each system may vary
- Transducer selection may not be aligned
- Bilateral/Unilateral settings may not be aligned
- Experience settings may not be aligned
- Quiet versus noise targets may not be aligned
- Proprietary Coupler Output for Flat Insertion Gain (CORFIG) Hawkins D, Cook J. Hearing aid software predictive gain values: how accurate are they? *Hear Jour.* 2003;56(7):26-34.

---

---

---

---

---

---

---

---

---

---

### DSL-Potential sources of inter-system variance- signals

- DSL provides SPL targets but varies them by type; speech, speech noise, pure tone
- Signals may be classified differently across systems
- DSL provides customized targets when the signals spectrum is provided to the algorithm
  - E.G. Otosuite feeds all available signal spectra to derive custom targets

---

---

---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

So why the difference in targets shown in the fitting software and the PMM system?

---

---

---

---

---

---

---

---

### Primary sources of difference between fitting and verification software

- Mismatch in understanding of what it means to fit to a prescriptive target versus selecting a prescription as the applied algorithm
- Mismatched signal types
- Mismatched versions of prescriptive target
- Mismatch in use of predicted versus measured REUG or RECD
- Mismatch in Binaural vs Monaural fitting selection
- Mismatch between Adult vs Pediatric or experience type

---

---

---

---

---

---

---

---

### What can I do about it

- Familiarize yourself with the prescription you are applying to your fittings
  - Keidser G, Dillon H, Flax M, Ching T and Brewer S (2011). The NAL-NL2 prescription procedure. *Audiology Research*, 1:e24:88-90.
  - Scollie, S., Seewald, R., Cornelisse, L., Moodie S., Bagatto, M., Lurnagaray, D., Beaulac, S., and Pumford, J. (2005). The Desired Sensation Level Multistage Input/Output Algorithm. *Trends in Amplification*, 9(4): pp. 159-197.
  - Johnson, E., & Dillon, H. (2011). A comparison of gain for adults from generic hearing aid prescriptive methods: Impacts on predicted loudness, frequency bandwidth, and speech intelligibility. *Journal of the American Academy of Audiology*, 22, 441-459.
- Check the settings in your PMM system and fitting software when you see a mismatch
- Consult reference manuals and clinical support from both

---

---

---

---

---

---

---

---



---

---

---

---

---

---

---

---

QUESTIONS???

---

---

---

---

---

---

---

---

What impact do some of the various settings have on the prescription?

---

---

---


---

---

---

---

---



---

---

---

---


---

---

---

---

### DSL



---

---

---

---

---

---

---

---

### NAL



---

---

---

---

---

---

---

---

So do I have to make selections in Fitting Details???

---

---

---

---

---

---

---

---



Why does the target sometimes fall at or below the clients hearing threshold levels?

---

---

---

---

---

---

---

---

Why does the target sometimes fall at or below the clients hearing threshold levels?

- Restoration of LTASS audibility of high frequencies is not always desirable
- In the case of NL2, the effective audibility factor it applies assumes that as hearing loss becomes more severe less information can be extracted from the speech signal
- Prescribed insertion gain may not aim to achieve audibility at higher frequencies

Keidser et al, 2011

---

---

---


---

---

---

---

---




---

---

---

---


---

---

---

---

Soft speech target and aided response for mild hearing loss



---

---

---

---

---

---

---

Thank you!!!

---

---

---

---

---

---

---